




### Vector Pyroelectric Detectors

Vector Pyroelectric Detectors, as with all Sciencetech detectors, have our traditional quality built in and offer futuristic features and performance. A variety of detectors are available including the unique high damage threshold HD and HDX models which can withstand energy densities up to 1.4 J/cm<sup>2</sup> and 12 J/cm<sup>2</sup> respectively. High resolution HR detectors can measure repetition rates up to 4 kHz and resolve energy readings to 3 nJ.



**SPECTRAL RESPONSE OF VECTOR PYROELECTRIC DETECTORS**

**FIGURE 1**

SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE.

The heart of the pyroelectric joulemeter, a sensitive detector of pulsed radiation, is the pyroelectric element. This element contains a temperature sensitive crystal which exhibits spontaneous electric polarization when exposed to the sensible heat generated by absorbed modulated radiation. The polarization is measured as a voltage at the electrodes which are attached to the crystal. This voltage is proportional to the laser energy level being absorbed. Since no power is required to drive the joulemeters, they can be directly connected to an oscilloscope as well as a Vector Synenergy Indicator. Futuristic Vector detectors are available with normal and slim profiles in highly absorbing, high rep rate, and high damage threshold models.

**VECTOR™ PYROELECTRIC DETECTOR SPECIFICATIONS**

Model No.	SP25	SPHF25	SPHD25	SP50	SPHF50	SPHD50	P25	PHF25	PHD25	PHDX25	PHDX25UV	PHDX50	PHDX50UV	P50	PHF50	PHD50
Aperture Size (dia.)	25.4 mm	25.4 mm	25.4 mm	50.8 mm	50.8 mm	50.8 mm	25.4 mm	25.4 mm	25.4 mm	7 mm	7mm	15 mm	15 mm	50.8 mm	50.8 mm	50.8 mm
Spectral Response	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1	See Fig. 1
Average Power (max.)	5 W*	5 W*	5 W*	10 W*	10 W*	10 W*	5 W*	5 W*	5 W*	5 W*	5 W*	10 W*	10 W*	10 W*	10 W*	10 W*
Minimum Energy	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**	7%**
Noise Equivalent Energy	4 µJ	4 µJ	4 µJ	16 µJ	16 µJ	16 µJ	4 µJ	4 µJ	4 µJ	4 µJ	4 µJ	16 µJ	16 µJ	16 µJ	16 µJ	16 µJ
Energy Density (max.)	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25	See Note 25
Accuracy	5%	5%	8%†	5%	5%	8%†	5%	5%	8%†	8%†	8%†	8%†	8%†	5%	5%	8%†
Output Sensitivity	8 V/J	8 V/J	2 V/J	2 V/J	2 V/J	2 V/J	8 V/J	8 V/J	2 V/J	2 V/J	2 V/J	2 V/J	2 V/J	2 V/J	2 V/J	2 V/J
Maximum Rep Rate	100 pps	400 pps	40 pps	50 pps	400 pps	20 pps	100 pps	400 pps	40 pps	40 pps	40 pps	40 pps	20 pps	50 pps	400 pps	20 pps
Maximum Pulse Duration	0.2 msec	0.045 msec	0.2 msec	0.4 msec	0.045 msec	0.4 msec	0.2 msec	0.045 msec	0.2 msec	0.2 msec	0.2 msec	0.4 msec	0.4 msec	0.4 msec	0.045 msec	0.4 msec
Dimensions (DxL) (cm.) (Width, SP Models)	5.8 x 5.8 x 1.4	5.8 x 5.8 x 1.4	5.8 x 5.8 x 1.4	7.6 x 7.6 x 1.5	7.6 x 7.6 x 1.5	7.6 x 7.6 x 1.5	6.12 x 5.77	6.12 x 5.77	6.12 x 5.77	6.12 x 9.88	6.12 x 8.89	8.76 x 8.88	8.76 x 9.88	8.76 x 5.77	8.76 x 5.77	8.76 x 5.77
Weight (kgs.)	0.14	0.14	0.14	0.18	0.18	0.18	0.41	0.41	0.41	0.5	0.5	0.77	0.77	0.68	0.68	0.68

\*Full illumination of sensor †Beam centered on absorber \*\*of selected range

**VECTOR HR PYROELECTRIC DETECTOR SPECIFICATIONS**


Model No.	PHF02	PHF05	PHF09	P05	P09
Aperture Size (dia.)	2 mm	5 mm	9 mm	5 mm	9 mm
Spectral Response	see fig. 1	see fig. 1	see fig. 1	see fig. 1	see fig. 1
Voltage Response (V/mJ)					
S, I	15	2.5	1	3.0	0.8
L	0.15	0.025	0.01		
Electrical Decay Time (msec) (RC Time Constant)				2.0	2.0
δ - pulse width ≤ 5 µsec	0.05	0.05	0.05		
I - pulse width ≤ 50 µsec	0.5	0.5	0.5		
L - pulse width ≤ 250 µsec	2.5	2.5	2.5		
Noise Equivalent Energy (nJ)				15	35
S, I	3.0	15	35		
L	150	750	3500		
Maximum Energy (µJ)				1.5	3.5
S, I	0.3	1.5	3.5		
L	15	75	350		
Maximum Rep Rate (pps)				400	200
S	4000	4000	4000		
I	400	400	400		
L	80	80	80		
Maximum Pulse Width (µsec) (For calibrated Response)				50	100
S	5	5	5		
I	50	50	50		
L	250	250	250		
Maximum Voltage Output (V)	4.5	4.5	4.5	4.5	4.5
Maximum Average Power (W)	1	2	2	2	2
Accuracy	7%	7%	7%	5%	5%
Maximum Energy Density	See Note 25				
Dimension (D x L) (cm.)	4.45 x 10.8				
Weight (kgs.)	0.21				
Indicator Compatibility	S310, S310D				

The P and SP models utilize a unique, durable, highly absorbing black coating which yields the flattest spectral response from the UV to the mid IR with reflection limited over the entire spectral region as indicated in Figure 1 (above) for the specifications. The high damage HD and HDX models are unique in pyroelectric measurement since they yield unsurpassed damage specifications. The HF high rep rate models have an energetic, partially absorbing, partially reflecting chromium coating. This coating is inherently fast and detects individual pulse energies at high repetition rates. The high resolution HR models are capable of low energy detection with noise equivalent energy levels down to 3 nJ on the most sensitive model. To accomplish accurate readings at such low energy levels, each HR detector contains a preamplifier powered by a 9 volt battery. This internal battery supplies noise-free power to the amplifying circuit. Inaccurate readings caused by low battery conditions are immediately and clearly recognized by a low battery warning LED light located on the top of the detector. All detectors come with a 1.27 cm (0.5") diameter x 8.89 cm (3.5") long post for mounting to an optical table and a BNC interconnect cable. A base (Model 301-019) for holding the detector/mounting post assembly upright is also available.

An optimum calibration of each detector is ensured by Sciencetech's highly trained and skilled service team. These experts use an electrically calibrated Sciencetech calorimeter/indicator system which has been NIST certified as the transfer standard. The energy of Sciencetech's YAG laser is accurately measured by the NIST certified calorimeter system and transferred to each Vector pyroelectric detector. This dependable calibration procedure generates a reliable certificate of calibration which accompanies each Vector joulemeter system and detector.

**ENERGY DENSITY FOR PYROELECTRIC DETECTORS**

Note 1	Max. J/cm <sup>2</sup> = 316 x (pulse width) <sup>1/2</sup> Maximum pulse width of the pyroelectric detector must be observed
Note 2	<b>HD Models</b> Max. J/cm <sup>2</sup> = 4500 x (pulse width) <sup>1/2</sup> to a maximum of 1.4 J/cm <sup>2</sup> Maximum pulse width of the pyroelectric detector must be observed
Note 3	<b>HDX Models</b> Max. J/cm <sup>2</sup> = 36,000 x (pulse width) <sup>1/2</sup> to a maximum of 12.6 J/cm <sup>2</sup> Maximum pulse width of the pyroelectric detector must be observed
Note 4	<b>HDXUV Models</b> Max. J/cm <sup>2</sup> = 36,000 x (pulse width) <sup>1/2</sup> to a maximum of 5.6 J/cm <sup>2</sup> Maximum pulse width of the pyroelectric detector must be observed



**MAXIMUM ENERGY DENSITY OF VECTOR PYROELECTRIC DETECTORS** **FIGURE 2**

Excessive energy or peak power density can damage any Vector pyroelectric detector. See Figure 2 (above) for safe operating limits. If your calculated specifications exceed those provided in Figure 2 and the specification table, optical attenuation is required. Scientech provides a damage test slide with each detector (except HD and HDX models) for preoperation damage testing.



**Features:**

- Joulemeter single pulse energy measurements at rep rates up to 4 kHz
- Advanced EMI/RFI protection
- High damage threshold models
- Slim profile models
- Easy connection to oscilloscope
- N.I.S.T. traceable calibration

Notes	
<b>Note 1</b>	AC2501, AC5001 30W/cm2@1064nm, 23W/cm2@532nm, 8.5W/cm2@355nm, 175mW/cm2@266nm
<b>Note 2</b>	ACX2501, ACX5001 Note 1 specs x 8 for 400nm to 1.2µm
<b>Note 3</b>	AC25UV, AC50UV, 384UV5, 388UV5, UC150UV 50W/cm2@355nm
<b>Note 4</b>	AC2504, AC5004 35W/cm2@1064nm
<b>Note 5</b>	AC2501, AC5001 100GW/cm2@1064nm, 78GW/cm2@532nm, 29GW/cm2@355nm, 580MW/cm2@266nm
<b>Note 6</b>	AC25UV, AC50UV, 384UV5, 388UV5, UC150UV Repetitive pulses: 101MW/cm2@355nm; Single pulses: 3.5GW/cm2@355nm
<b>Note 7</b>	AC2504, AC5004 125GW/cm2@1064nm
<b>Note 8</b>	AC2500, AC5000, 360401, 360801, UC150 Max J/cm2 = 1,000 x (pulse width)1/2 to a max of 200J/cm2
<b>Note 9</b>	AC250FX, AC50FX Max J/cm2 = 4,950 x (pulse width)1/2 to a max of 12.3J/cm2
<b>Note 10</b>	ACX25FX, ACX50FX Max J/cm2 = 39,600 x (pulse width)1/2 to a max of 36.9J/cm2
<b>Note 11</b>	AC2501, AC5001 Repetitive pulses: 4.1J/cm2@1064nm, 3.2J/cm2@532nm, 1.2J/cm2@355nm, 24mJ/cm2@266nm Single pulses: 8J/cm2@1064nm, 6.2J/cm2@532nm, 2.3J/cm2@355nm, 46mJ/cm2@266nm
<b>Note 12</b>	ACX2501, ACX5001 Note 11 specs x 8 for 400nm to 1.2µm
<b>Note 13</b>	AC25UV, AC50UV, 384UV5, 388UV5, UC150UV Repetitive pulses: 1.1J/cm2@355nm; Single pulses: 40J/cm2@355nm
<b>Note 14</b>	AC25004, AC5004 Repetitive pulses: 4.8J/cm2@1064nm; Single pulses: 10J/cm2@1064nm
<b>Note 15</b>	AC25FX, AC50FX Max p.p.d.: 70MW/cm2@1064nm pulse
<b>Note 16</b>	ACX25FX, ACX50FX Max p.p.d.: 560MW/cm2@1064nm pulse
<b>Note 17</b>	380401 27W/cm2@1064nm, 21W/cm2@532nm, 7.7W/cm2@355nm, 158mW/cm2@266nm
<b>Note 18</b>	380801 13.5W/cm2@1064nm, 10.5W/cm2@532nm, 3.85W/cm2@355nm, 79mW/cm2@266nm
<b>Note 19</b>	380401 90GW/cm2@1064nm, 71GW/cm2@532nm, 27GW/cm2@355nm, 530MW/cm2@266nm
<b>Note 20</b>	380801 45GW/cm2@1064nm, 35.5GW/cm2@532nm, 13.5GW/cm2@355nm, 265MW/cm2@266nm
<b>Note 21</b>	380401 Repetitive pulses: 3.7J/cm2@1064nm, 2.9J/cm2@532nm, 1J/cm2@355nm, 20mJ/cm2@266nm Single pulses: 7J/cm2@1064nm, 5.6J/cm2@532nm, 2.1J/cm2@355nm, 41mJ/cm2@266nm
<b>Note 22</b>	380801 Repetitive pulses: 1.85J/cm2@1064nm, 1.45J/cm2@532nm, 0.5J/cm2@355nm, 10mJ/cm2@266nm Single pulses: 3.5J/cm2@1064nm, 2.8J/cm2@532nm, 1.05J/cm2@355nm, 20.5mJ/cm2@266nm
<b>Note 23</b>	AC25HD, AC50HD, UC150HD, UC150HD40 Max J/cm2 = 4,500 x (pulse width)1/2 to a max of 14J/cm2
<b>Note 24</b>	ACX25HD, ACX50HD Max J/cm2 = 36,000 x (pulse width)1/2 to a max of 42.5J/cm2
<b>Note 25</b>	SP25, SPHF25, SP50, SPHF50, P25, PHF25, P50, PHF50, PHF02, PHF05, PHF09, P05, P09 Max J/cm2 = 316 x (pulse width)1/2
<b>Note 26</b>	SPHD25, SPHD50, PHD25, PHD50 Max J/cm2 = 4,500 x (pulse width)1/2 to a max of 1.4J/cm2
<b>Note 27</b>	PHDX25, PHDX50 Max J/cm2 = 36,000 x (pulse width)1/2 to a max of 12.6J/cm2
<b>Note 28</b>	PHDX25UV, PHDX50UV Max J/cm2 = 18,000 x (pulse width)1/2 to a max of 5.6J/cm2