

HYPHO 35W



2.1 μm radiation of the HyPHo system is strongly absorbed by water and tissue. Minimal penetration, less than 0.5mm, ensures safe tissue ablation. The laser operates only in pulsed mode. Four operative modes are available: Ablation, Coagulation, Fragmentation mode, Dusting mode. Tens of preset protocols are included in the software.

- Unique performances
- Unparalleled versatility
- Four operative modes
- Dedicated smart optical fibers
- Remote diagnostic
- Clinical data sharing

APPLICATIONS

General and percutaneous lithotripsy (also Pediatric lithotripsy):

- Renal calculi
- Ureteral calculi
- Bladder calculi

Urology: Open and endoscopic surgery

- Urethral strictures
- Ureteral strictures
- Bladder neck incisions
- Ablation and resection of bladder, urethral, ureteral and Kidney tumors
- Transurethral incision of the prostate
- Condylomas
- Lesions of external genitalia

Gastroenterology

Arthroscopy

Pulmonary

Gynecology

General Surgery

ENT



Request brochure

